

# MagniLink Go!



# MagniLink Go!

Imagine getting up in the morning, having a cup of coffee and reading the daily mail. For most people this is just a normal daily routine. For people with visual impairment reading the daily mail would involve a video magnifier or a reading machine. Now a faster and portable solution is here, making information accessible anywhere and whenever necessary - MagniLink Go!

## Become more independent!

MagniLink Go! is a new portable solution from LVI for your cell phone, allowing you to read whatever you want, wherever you might be. It turns your cell phone into a combined reading machine and digital magnifying glass, but it is really a software for your cell phone and accessibility on the go is now a reality.

## A single key press - and you are ready to Go!

MagniLink Go! gives immediate access to any printed material, such as a letter. It uses the cell phone's camera to either magnify a text using the live camera option or capture an image for speech output.

Choose a document, place the cell phone above the document and simply press a single button to begin the OCR processing. The speech will start automatically and users can listen to the content of the document aloud in a high quality voice, and then choose the preferred speech volume and speed. Users who can

see the screen can enlarge, choose color combination and choose between word or column scroll. Each function only requires a single key press, so there's no need to remember complicated short keys or go through difficult applications.

## Constantly on the go - use MagniLink Go!

Just a simple press - and you are ready to go and access all information wherever you might be. At the airport to receive gate information, in a restaurant to choose from a menu, or maybe to access the information in the daily mail while drinking your coffee in the morning.

It is here, and it is MagniLink Go!



## Features & benefits

- MagniLink Go! combines a magnifying glass and reading device - in your cell phone. This is useful for low vision users, locate what information you want with the live camera function and select it before you capture the picture.
- Capture a picture of a newspaper and then send the text to your friends using SMS, MMS or E-mail.
- The MagniLink Go! can read any document regardless of the camera's angle in relation to the paper. This provides an even faster process; just direct the camera to the text and the software will automatically adjust the angle.
- Any image file can be saved on the phone and then OCR processed. For instance receive a MMS of a plane ticket with gate number, and then turn the information into speech.

## Further possibilities

- Combine MagniLink Go! with MobileSpeak or Talks for an all-in-one mobile solution. With the use of a screen reader, the user can access the other functions of the phone device including making and receiving phone calls, personal information management such as contacts, calendar, etc.
- Select different accessories from Nokia for your cell phone, such as headphones, Bluetooth external speakers or Bluetooth external keyboards.



accessible  
anywhere and  
whenever  
necessary -  
MagniLink Go!



# TECHNICAL SPECIFICATIONS

Model	MLGO
Supported phones	Nokia N82, N86 6220 classic (read more at <a href="http://www.lvi.se">www.lvi.se</a> )
System requirements	Operating system: Symbian 9.3 S60. 3rd Ed, FP1 and FP2 Camera: Minimum 5 megapixel with autofocus
Voices	Acapela: High Quality Text To Speech (TTS) (Choose preferred voices when purchasing the software) Nokia: Built-in voices from Nokia
Application languages	English, Spanish, Danish, Dutch, Finish, French, German, Norwegian, and Swedish
Magnification	1 - 7 x
Color modes	Natural colors, 6 different artificial colors (positive/negative)
Updates	Free software updates within version
Accessories	2 polarizing filters and a cell phone strap are included

Northern California Dealer:



Showroom Location  
1400 Terra Bella Avenue, Suite I  
Mountain View, CA 94043  
Tel: 650-961-6541  
Fax: 650-968-4740  
[www.freedomvision.com](http://www.freedomvision.com)

A Product Of:

